## **BondSaver**

## **BondSaver**



**BondSaver** has been developed to provide 5 significant benefits:

- Desensitizes Instantly
- Prevents MMP Generation and Bond Degradation
- Disinfects with Potent Antimicrobial Agents
- Stiffens and Fixes Collagen
- •Provides Fluoride at the Restorative-Tooth Interface
- •Easy to apply, no drying, no rinsing, no scrubbing.... Simply apply and move on in your protocol.

Simply put, BondSaver optimizes the Restorative-Tooth interface.

## **BONDSAVER versus DESENSIMAX**

BONDSAVER has been formulated to be used with self etch adhesives, glass ionomers, resin modified glass ionomers.

DESENSIMAX has been formulated to be used with total etch and selective etch adhesive systems along with the treatment of sensitivity resulting from exposed dentin.

