## RE-GEN SELF-ETCH ADHESIVE



## **RE-GEN™ Self-Etch Adhesive**



RE-GEN™ Self-Etch is a 3-bottle system that consists of a conditioner, primer and bonding resin. RE-GEN™ combines clinically proven adhesives with the power of Bioglass, a material verified to stimulate bone cells towards a path of regeneration and self-repair. Bioactive material that forms apatite on the surface would have several benefits, including the closure of any gaps forming at the resin-dentin interface, and potentially improved long-term bond strengths, eliminating bond degradation.

- •The release of calcium, phosphate and fluoride to promote hydroxyapatite formation
- •Bond Strength: 55MPa
- •Efficient application protocol takes only 35 seconds
  - Only one protocol for all procedures
- Can be used for both direct and indirect restorations
- No-rinse system eliminates any post-operative patient sensitivity
- •No technique sensitivity promotes ease-of-use and consistent results
- •Contains No Bis-GMA, No Bisphenol A, No BPA derivatives

